




Distortion in STEM Spectrum imaging (SI) ?




It works for 2 types of distortion in STEM SI:


(a) linear distortion




(b) nonlinear distortion




How to install and use?

Copy the gtk file into the DM plugin folder. (for the nonlinear a HREM mouse plugin is required )




The scripts will appear under [STEM-SI Warp] menu:



(a) For linear distortion:

1. Select script “**linear warp**”
2. Input the angle α [Positive: Warp to left; Negative Warp to right]
4. Copy the ROI, and crop ROI



This Script works for 2D - 4D images!






(b) For nonlinear distortion:

1. Use HREM mouse tool to select the atomic columns that you want to align to a straight line



2. Select script "nonlinear warp -2D", Image 0: select the marked ADF image, Image 1 and 2: the element maps




If you want to warp multi-dimensional data, select 3D or 4D accordingly !


4. Repeat on other elements maps

That's it !

3. Input the pixel size for atomic column locating



You will get warped ADF image and element map



Use copy ROI and crop ROI script to Crop the area of interest.

